

## **Key Results**

- EUV photon energy in the range of **30-120eV** pumped with 785nm wavelength.
- Bright harmonics with narrow linewidths pumped at a fixed 1kHz repetition rate at a customer laboratory. Up to 10kHz repetition rate can also be pre-selected.
- High flux EUV photon generation observed using a minimum pump pulse energies of 0.5mJ with Argon and of 2mJ with Helium.
- High EUV photon flux at source is measured and calibrated, reaching up to 8.0x 10<sup>11</sup> photons/sec full spectrum and 1.7x 10<sup>11</sup> photons/sec per harmonic at the 27<sup>th</sup> harmonic (43eV in photon energy).
- Driven with <40fs pulses, the EUV pulses are suitable for a wide range of applications, including but not limited to, spin dynamics, attosecond nonlinear physics, trARPES, and more.
- The XUUS-5 and Arterium beamline installation were done within days and all the HHG performance measurements were done in-situ. The integration process is quick and easy.
- The resulting source configuration is simple, versatile, commercially available and readilyaccessible to customers who have an equivalent Ti:sapphire pump source.

## **HHG Performance**

|                           | Ar                    | He                    |
|---------------------------|-----------------------|-----------------------|
| Rep Rate (kHz)            | 1                     | 1                     |
| Pulse Energy (mJ)         | 1.0                   | 4.0                   |
| Harmonic of 785nm         | 27                    | 61                    |
| Photon Energy (ev)        | 43                    | 96                    |
| Flux At Source (phot/sec) | 8.0x 10 <sup>11</sup> | 2.4x 10 <sup>10</sup> |

Representative flux measurements with Si diodes for all harmonics and calculated to source based on 50% rejector reflectivity for 30nm and 80% for 13nm, and 2-3% filter transmission. Single harmonics range from 10-40% of total flux.

## **HHG** spectra





